salling

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Vi tager forbehold for trykfejl.

INTRODUCTION

Thank you for choosing a product from Salling. You have acquired a quality product which meets all applicable performance and safety standards in the EU. In order to ensure proper treatment and a long product life, we recommend that you comply with the instructions below.

Please read the manual and especially the safety instructions carefully before using the product. Keep this manual and provide it for referral along with the product when others are using it.

GENERAL SAFETY INSTRUCTIONS

- Improper use of the product can cause personal injury and damage.
- Only use the product for what it is intended. The manufacturer is not liable for damages arising from improper use or handling.
- Always place the product on a dry, flat and stable surface well away from table edge or similar.
- The product is intended for household use only. The product is not for professional use.
- The product engine/electronics, cord or plug cannot be submerged in water, and no water may get in contact with the engine/electronics.
- Never touch the product, cord or plug with wet or damp hands.
- Do not expose the product to direct sunlight, high temperatures, humidity, dust or corrosive substances.
- Never leave the product when it is turned on.
- When the product is in use, it should be kept under constant supervision. Children should always be supervised when using the product to ensure that they do not play with the product. The product is not a toy.
- The product may have sharp parts and must, therefore, be kept out of reach of children.

- Only use the accessories that come with the product or that are recommended by the manufacturer.
- The product must not be used by people with reduced sensitivity, physical or mental disability, or people who are not able to operate the product unless they are supervised or instructed by a person responsible for their safety.

CORD AND PLUG

- Do not let the cord hang over the edge.
- Keep the cord and the product away from heat sources, hot objects and open flames.
- Unplug the power cord before cleaning and if the product is not in use. Do not pull the cord when removing the plug from the socket, rather grasp the plug itself.
- Make sure that the cord is fully extended.
- Check regularly that the cord and the plug are not damaged and do not use the product if this is the case, or if it has been dropped on the floor, dropped in water or damaged in another way.
- If the product or plug is damaged, the product must be inspected and, if necessary, repaired by an authorised repairer. Never try to repair the product yourself. Please, contact the store to see if there is a repair that falls under warranty. Tampering of the product will void the warranty.

RULES OF COMPLAINTS THE WARRANTY DOES NOT APPLY:

- If the above is not observed.
- If there has been unauthorised intervention of the product.
- If the product has been mishandled, subjected to rough treatment or some other form of damage.

PRODUCT DISPOSAL

Used product should not be disposed as household waste. According to legislation, the product should be disposed of properly so that the materials can be recycled in an environmentally

friendly manner. Further information is available from the responsible local authority or local disposal company. The product is equipped with the following logo:

Danger: In order to reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions.

For adults only

It can only be used on dry, clean skin surfaces of the body, and cannot be used through clothes. Press and move it gently on the skin. The duration on each part is about 60 seconds. Use this device only on soft tissues of the body without causing pain or discomfort. Do not use it on head or any hard or boned part of the body.

Only use the vibrating head that can produce the best results.

Bruise may be caused whether under moderate pressure or normal pressure. Please always check the area that receives physiotherapy.

Stop using it immediately once there is any sign of pain or discomfort.

Keep fingers, hair, or other parts of the body away from the shaft and back of the vibrating head, because they might be pinched.

Do not place any object in the vent of the device.

Do not immerse it in water or allow water to enter the vent of the device.

Do not drop the device or use it improperly.

Charge only with DC 5V input (adapter not included)..

Carefully check the device and battery before each use.

Do not change or modify the device in any way.

Never operate or charge the device in an unattended way.

WARNINGS AND REMINDERS

Do not immerse it into water. Keep it away from liquids or heat sources. Keep the vent free from dust and scraps.

Do not remove the screws or try to disassemble the device.

Do not continuously use it for more than half an hour. Please stop the device for 15 minutes and then use it again.

Remove the USB cable (this product comes with a USB cable only, adapter not included) from the massager after charging or before use.

It is only for adults. If you have injuries, do not use this device. Please consult with your doctor before using this device.

CONFIGURATION AND SPECIFICATION

Configuration

- 1. Power on/OFF button
 - 2. Power level indicator lamp LED (blue)
 - 3. Handle/fix battery
 - 4. Removable, replaceable vibration head
- When removing the vibration head, just exert your strength to pull it out from the main appliance.
- 5. Type-C cable
 - 6-level adjustment meets different needs

Level indicator lamp

- Level 1 - 2500 rpm
- Level 2 - 2700 rpm
- Level 3 - 2900 rpm
- Level 4 - 3100 rpm
- Level 5 - 3300 rpm
- Level 6 - 3500 rpm
- Level adjustment switch
- Charging hole

U head/for vertebrae

Small round head/ for joints

Ball head/for large muscle groups

Small flat head/for all parts of the body

MAINTENANCE, CLEANING, AND STORAGE:

Clean the device with a slightly wet towel, and then wipe it dry with a soft cloth for storage or put it into the supplied carton when you go traveling.

CHARGING

Before using it for the first time, charge it fully for three hours. Connect to DC5V input (adapter not included) for charging.

When the battery is connected, the LED on the device flashes and displays the battery level, indicating that the battery is being charged.

When all LEDs stop flashing and are on, it indicates that the battery is fully charged.

It is not recommended to use up the battery to the lowest blue LED level.

The average operating time of the device is half an hour, depending on the speed and the applied pressure during the use.

TECHNICAL DATA

Article No. 10958217

Input: DC 5V, 1A

25 Watt

2 x 800 mAh Li-ion Battery

PRODUCER

Epiq ApS

Krøyer Kielbergs Vej 3, 3. th

8660 Skanderborg, Denmark

www.epiq.dk

We cannot be held responsible for any printing errors.