# Model 74136

Brugsanvisning

Bruksanvisning

Bruksanvisning

Käyttöohje

Instruction manual

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Kasutusjuhend

Manual de instrucciones

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Manuel d'instructions

DK

NO

SE

FI

GB

DE

PL

ET

ES

IT

NL

FR

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# NECK AND SHOULDER MASSAGER

## Introduction

To get the most out of your new neck and shoulder massager, please read these instructions before use. Please also save the instructions in case you need to refer to them at a later date.

## Technical data

Voltage/frequency: 12 V DC/2.0 A

Power rating: 25 W

With 8 rotating Shiatsu massage heads

With heat and knead function, rotation clockwise and anticlockwise, and automatic shutdown after 15 minutes

Includes adapters for 230 V and 12 V DC (cigarette lighter socket)

Class III product: Designed to be supplied from a separated extra low voltage supply (SELV). The voltage from a SELV supply is low enough that a person can normally come into contact with it without the risk of electric shock.

Do not use this appliance in a bathtub, shower, or water-filled bowl.

## Special safety instructions

This appliance must not be used by persons, including children, with a physical, sensory or mental disability or anyone without experience or knowledge unless they are supervised or have been instructed how to use the appliance.

Children should be supervised to ensure that they do not play with the appliance.

If you are pregnant, have a pacemaker or have diabetes or any other health issues, you should ask your doctor for advice before using the massager.

If you have persistent muscle or joint pain that does not disappear after using the massager, you should ask your doctor for advice.

The massager's heat function must not be used by people who are particularly sensitive to heat, or who find it difficult to feel heat.

Do not use the massager on areas of skin with open wounds, or if the skin is discoloured, burned, inflamed or affected by eczema.

The massager should not be used for more than 15 minutes at a time on the same muscles, and for no more than 30 minutes at a time in total.

Do not use the massager while you are lying in bed or on a sofa.

Check regularly that the massager, cord and plug are in perfect condition. If this is not the case, the massager must not be used.

Do not attempt to repair the massager yourself.

The massager, cord or plug must not be exposed to heat, direct sunlight, moisture, sharp edges and the like.

Always turn the massager off and unplug it from the socket when not in use or when cleaning it. Grasp the plug to pull it out, not the cord.

The massager is solely intended for domestic use and is not suitable for commercial use. Always use the massager in accordance with the instructions.

## Main components

1. Handle
2. Massage heads
3. Connector for power supply
4. Button for heat function
5. On/off button
6. Button for direction of rotation
7. 12 V DC adapter with cigarette lighter connector
8. 230 V adapter

GB

## Use

Connect the required adapter (7, 8) to the socket (3) on the massager.

Plug the adapter (7) into a 12 V cigarette lighter socket in a car or the adapter (8) into an ordinary residential mains socket.

Place the massager around the neck and place your hands in the handles (1).

Switch the massager on using the on/off button (5). The built-in massage heads start to rotate against your neck and shoulders.

Press the button (4) to switch the heat function on or off.

Press the button (6) to select the direction of rotation.

Pull on the handles to control the massage pressure.

Switch the massager off by pressing the on/off button once you have finished using it. If you have not already switched the massager off, it will turn off automatically after 15 minutes of massage.

Unplug the adapter from the socket on the massager.

## Cleaning and maintenance

Wipe the massager with a dry cloth.

The massager must not be washed or exposed to water.

Do not use corrosive or abrasive cleaning agents.

Store the massager in a dry, dust-free location out of direct sunlight.

## Service centre

**Note: Please quote the product model number in connection with all inquiries.**

The model number is shown on the front of this manual and on the product rating plate.

For:

- Complaints
- Replacement parts
- Returns
- Guarantee issues
- [www.schou.com](http://www.schou.com)

## Environmental information

Electrical and electronic equipment (EEE) contains materials, components and substances that may be hazardous and harmful to human health and the environment

if waste electrical and electronic equipment (WEEE) is not disposed of correctly. Products marked with a crossed-out wheeled bin are electrical and electronic equipment. The crossed-out wheeled bin indicates that waste electrical and electronic equipment must not be disposed of with unsorted household waste, but must be collected separately.

Manufactured in P.R.C.

Manufacturer:  
Schou Company A/S  
Nordager 31  
DK-6000 Kolding

All rights reserved. The content of this manual may not be reproduced, either in full or in part, in any way by electronic or mechanical means, e.g. photocopying or publication, translated or saved in an information storage and retrieval system without written permission from Schou Company A/S.

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