TUNTURI®

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GB	User manual	06 - 10
DE	Benutzerhandbuch	11 - 15
FR	Manuel de l'utilisateur	16 - 20
NL	Gebruikershandleiding	21 - 25
IT	Manuale d'uso	26 - 30
ES	Manual del usuario	31 - 35
SV	Bruksanvisning	36 - 40
SU	Käyttöohje	41 - 45

Attention, Achtung, Attention, Attentie, Attenzione, Atención, Giv akt, Huomio

- Please read this user manual carefully prior to using this product.
- Lesen Sie sich vor dem Gebrauch des Produkts diese Bedienungsanleitung sorgfältig durch.
- Veuillez lire attentivement le présent manuel de l'utilisateur avant d'utiliser ce produit.
- Lees deze gebruikershandleiding zorgvuldig voordat u dit product gaat gebruiken.
- Leggere attentamente questo manuale dell'utente prima di utilizzare il prodotto.
- Antes de utilizar este producto, lea atentamente este manual de usuario.
- Läs den här användarhandboken noggrant före användning av den här produkten.
- Lue tämä käyttöohje huolellisesti ennen tämän tuotteen käyttöä.

A

B

C

D-1

D-2

D-3

D-4

Index

- Safety warnings..... 6**
- Control panel (Fig. B) 7**
- Box content (Fig. C)..... 7**
- Usage instructions 7**
- Charging..... 7
- Operating 7**
- Power on / off 7
- Speed levels 8
- Maintenance, cleaning & storage 8**
- Troubleshooting 8
- Defects & malfunctions..... 8
- Specifications..... 9**
- Massage head user instructions 9**
- Warranty..... 9**
- Declaration of the manufacturer 10**
- Disclaimer 10**

Massage Gun

Safety warnings

WARNING

- Read the safety warnings and the instructions. Failure to follow the safety warnings and the instructions can cause personal injury or damage to the equipment. Keep the safety warnings and the instructions for future reference.

DANGER

- To reduce risks of electric shock, fire and personal injury, this product must be used in accordance with the following instructions.

- The equipment is suitable for domestic use only. The equipment is not suitable for commercial use.
- The use of this equipment by children or people with a physical, sensory, mental or motorial disability, or lack of experience and knowledge can give cause to hazards. People responsible for their safety must give explicit instructions or supervise the use of the equipment.
- If you experience nausea, dizziness or other abnormal symptoms, immediately stop your workout and consult a physician.
- The equipment is suitable for indoor use only. The equipment is not suitable for outdoor use.
- Only use the equipment in environments with ambient temperatures between 10°C and 35°C. Only store the equipment in environments with ambient temperatures between 5°C and 45°C.
- Do not use or store the equipment in humid surroundings. The air humidity must never be more than 80%.
- Only use the equipment for its intended purpose. Do not use the equipment for other purposes than described in the manual.
- Do not use the equipment if any part is damaged or defective. If a part is damaged or defective, contact your dealer.
- Keep your hair away from the moving parts.
- Keep jewellery and other objects away from the moving parts.
- Keep fingers, hair or any other body part away from the shaft and backside of the massage head as pinching may occur.
- Make sure that only one person uses the equipment at a time.
- Do not open equipment without consulting your dealer.
- Use only on a dry clean surface of the body by lightly pressing and moving across the body for approximately 60 seconds per region.
- Use the Massage Gun only on the body's soft tissue as desired without producing pain or discomfort.
- Do not use on head or any hard or bony area of the body.
- Use only the massage heads that provide the best desired effect.
- Bruising can occur regardless of control setting or pressure applied. Check treated areas frequently and immediately stop at first sign of pain or discomfort
- Do not place any objects into the ventilation holes of the Massage Gun.

- Do not immerse in water or allow water to enter the ventilation holes of the Massage Gun.
- Do not drop or misuse the Massage Gun.
- Only recharge the massage gun with a 5 volt USB charger in combination with an USB-C cable.
- Carefully examine device and battery before each use.
- Do not tamper or alter Massage Gun in any way.
- Never leave the Massage Gun operating or charging unattended.
- Keep ventilation ports free from dust and debris.
- Do not operate continuously for more than one hour. Allow device to rest for 30 minutes before re-using.
- Unplug the charger from device battery after charging or prior to use.
- Do not use if injured. Consult your doctor before using this product.

Control panel (Fig. B)

Your Tunturi Massage Gun Mini is equipped with a control panel on the bottom side of the hand grip.

1. Speed level indicator LED's
2. Control Button.
3. USB-C Charging port

Box content (Fig. C)

Your Tunturi Massage Gun Mini comes with the following items in the carton.

1. Storage case
2. Massage gun mini
3. 4 massage heads
4. USB charging cable
5. User manual

Usage instructions

Charging

- Fully charge battery for at least 3 hours before first use.
- To charge, connect the USB cable type-C end of the supplied cable to the USB charging port on the massage gun handgrip's bottom end. Connect the USB cable Type-A end of the supplied cable to a USB charger and plug the charger into a power socket.
- LED ring on handgrip lights pulsate to show battery level and indicate active charging.
- Full charge is indicated when the LED ring light stops flashing and remains fully illuminated.
- Battery may be recharged at any time and at any battery level.
- It is not recommended to fully deplete the battery to the red LED level.

NOTE

- When the battery power is less than 20%, the light ring at the bottom of the machine shows lights red.

Operating

The Tunturi Massage Gun is a portable rechargeable massage device with large torque and long adjustment. It can effectively reduce the pain caused by the accumulation of lactic acid after intense exercise, and relieve the discomfort of muscles in the sedentary and long-lying population.

Power on / off

- To power on the device press and hold the control button on bottom panel for 3 seconds when being in off mode.
- Green LED ring will light up and illumines constantly. (Standby mode)
- To power off the device press and hold the control button on bottom panel for 3 seconds.

NOTE

- The massage gun cannot be turned on when being in charging mode.
- The massage gun will switch into "Charging mode" when connected to a power source.

Speed levels

Use the control button to toggle between the different speed levels.

- Press control button shortly to change speed level.
- The Speed level indicator LED's shows the operation speed of the massage gun.

Toggle sequence is :

Speed 1	1 LED lights up	1400 rpm
Speed 2	2 LED lights up	1680 rpm
Speed 3	3 LED lights up	1960 rpm
Speed 4	4 LED lights up	2240 rpm
Speed 5	5 LED lights up	2520 rpm
Speed 6	6 LED lights up	2800 rpm
Standby	0 LEDs lights up	0 rpm

NOTE

- By pressing and holding the control button in operation mode, you will switch the massage gun in Off mode.

Maintenance, cleaning & storage

- Use a slightly damp towel to clean device and a dry with a soft cloth.
- Avoid moist to enter the device's housing at all times.
- For storage or travel store in supplied storage box.

NOTE

- When not using the massage gun for a longer period of time, we recommend you to charge the device at least once a month. This is to keep the battery in optimal condition.
- An exhausted battery will shorten it's lifetime or cause performance failure.

Troubleshooting

The massage gun does not have an error code signalling system. When there is too much stress applied to the massage gun when being in use, the massage gun will shut off.

Try to resolve by:

- When the massage gun is hot, give it some time to cool off.
- Press the control button and hold for 3 seconds to turn in Stand By mode again, and operate it as normal.

When a abnormal sounds comes out of the massage gun during usage, stop using the massage gun to avoid damage.

Try to solve by:

- Reinstalling the massage head to make sure the fixture is correct.
- Applying less load on the massage gun when using.
- When the massage gun will not turn on, try to charge it for at least 4 hours to make sure the battery is charged.

Defects & malfunctions

Despite continuous quality control, defects and malfunctions caused by individual parts can occur in the equipment. In most cases, it will be sufficient to replace the defective part.

- If the equipment does not function properly, immediately contact the dealer.
- Provide the model number and the serial number of the equipment to the dealer. State the nature of the problem, conditions of use and purchase date.

Specifications

Model	Message Gun mini
Net Weight	460 gr.
Battery Capacity	1500 M.Ah.
Power	12 W.
Dimensions	100x 47x 155 mm.
Charger power output	5 V-DC./ 2.0 A
Amplitude	6 mm
Charge time	3 hours
Full charge battery use time at:	
Level 6 speed	± 2.0 hours
Level 3 speed	± 2.5 hours
Level 1 speed	± 3.0 hours

Massage head user instructions

With unit switched in OFF mode, fully insert massage head by pressing in firmly.
With unit switched in OFF mode, remove the massage head by pulling it out. To ease the release you can rotate slightly and pull at the same time.

TIP

- When remove the massage head it is best to pull and turn the massage head, as if it is a threaded connection, to make it come loose more easily.
- When place the massage head it is best to push and turn the massage head, as if it is a threaded connection, to make secure the fixation more easily.

WARNING

- Do not press/ pull the massage head in/ out at an angle to avoid damages.

Round shape (fig D-1)

Round spherical design produces a power output of ropagating and diffusing type between receiving and pushing. Strength is even and soft. It is suitable for massage and relaxation of waist, abdomen, buttocks and legs. As well as muscle relaxation before and after exercise, recommend light weight fitness enthusiasts like dancing, yoga, aerobics and so on.

Mushroom shape (fig D-2)

Small and brief design, on the basis of the plane, a small arc which fits the body modelling mechanism is added and make it feel better. Suitable for the inside and outside muscles of the abdomen, chest and legs of the human body.

U Shape (fig D-3)

The force path between the two points is even, the stimulation is in place and powerful, which can effectively carry out some simple acupoint massage. Reasonable U-shaped width size perfectly avoids the human spine and cervical vertebra, and achieves a better sense of experience. Suitable for shoulder, cervical spine and bilateral muscles. (Can not directly hit the cervical spine and spine, it will cause physical injury).

Conical shape (fig D-4)

Simple shape, full of impact experience, direct pressure to the depths. Suitable for meridians, soles of feet, palms, Achilles tendons, rhomboid muscles of back, etc.

Warranty

Tunturi purchaser's warranty

Warranty terms

The consumer is entitled to the applicable legal rights stated in the national legislation concerning the commerce of consumer goods. This warranty does not restrict these rights. The Purchaser's Warranty is only valid if the item is used in an environment approved by Tunturi New Fitness BV for that particular equipment. The product-specific approved environment is stated in the Owner's Manual provided with your equipment.

Warranty conditions

The warranty conditions commence from the date of purchase. Warranty conditions may vary per country, so please consult your local dealer for the warranty conditions.

Warranty coverage

Under no circumstances shall Tunturi New Fitness BV or the Tunturi Distributor be liable by virtue of this warranty or otherwise for any special, indirect, secondary or consequential damages of any nature arising out of any use or inability to use this equipment.

This warranty excludes.

- Damage caused by accident, abuse, mishandling, or transport;
- Units not used in accordance with instructions;
- Damage exceeding the cost of the product;
- Deterioration of the delivered product resulting from abnormal storage and / or safeguarding conditions on the client's premises;
- Failure to provide the dated proof of purchase.

Failure to follow the instructions given in the Owner's Manual will invalidate the product's warranty.

Declaration of the manufacturer

Tunturi New Fitness BV declares that the product is in conformity with the following standards and directives: EN 957 (HC), 2014/30/EU. The product therefore carries the CE label.

03-2021

Tunturi New Fitness BV
Purmerweg 1
1311 XE Almere
The Netherlands

Disclaimer

© 2021 Tunturi New Fitness BV
All rights reserved.

- The product and the manual are subject to change.
- Specifications can be changed without further notice.
- Check our website for the latest user manual version.

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